



We have released a new version! The release is ready for download on our FTP-site. The web site is updated with new download links, and we posted an announcement.

Note that you will need to have an up-to-date license to install this update.

Changelog updates

- Manual and silent installations with MSI update/setup
- Support VDDS linking (images stages 1-5 implemented)
- Improved support child tooth positions
- Device added: Dexis Titanium sensor
- Device added: Kavo IXS sensor
- Image acquisition now possible for child tooth positions
- Support for Azure cloud MSSQL databases
- Improved import for KAVO/Palodex NCI interface
- Device added: Cruxell Cruxcan
- Image view added: 4x3 grid
- Integration with 3Shape Dental Desktop
- Integration with Boneprox portal
- Device added to remote Xray: FireCR
- Support for Firebird database

Changelog minor fixes

- Fixed issue minor issues.